

CHANGZHOU MEDITECH TECHNOLOGY CO., LTD

ADD:9th floor, INNO LASER Building, Science & Education Town, Wujin District, Changzhou, Jiangsu, China

WhatsApp: +86-18112515727

FAX: 86-519-85855955

Email: song@orthopedic-china.com

MT Calcaneal H-Locking Plate

<i>Products</i>	<i>Picture</i>	<i>REF</i>	<i>Specification</i>	<i>Thickness</i>	<i>Width</i>	<i>Length</i>
MT Calcaneal H -Locking Plate(Use 3. 5 Locking Screw)		5100-9701	6mm Height	1.8	22	25
		5100-9702	8mm Height	1.8	22	25
		5100-9703	10mm Height	1.8	22	25

CHANGZHOU MEDITECH TECHNOLOGY CO., LTD

ADD:9th floor, INNO LASER Building, Science & Education Town, Wujin District, Changzhou, Jiangsu, China

WhatsApp: +86-18112515727

FAX: 86-519-85855955

Email: song@orthopedic-china.com
