

CHANGZHOU MEDITECH TECHNOLOGY CO., LTD

ADD:9th floor, INNO LASER Building, Science & Education Town, Wujin District, Changzhou, Jiangsu, China

WhatsApp: +86-18112515727

FAX: 86-519-85855955

Email: song@orthopedic-china.com

Straight Type With Middle Pins External Fixator

<i>Products</i>	<i>Picture</i>	<i>REF</i>	<i>Specification</i>	
Straight Type with Middle Pins External Fixator		9100-0401	L 285mm	6*130mm Pin 6 PCS 6mm Hex Wrench 6mm Screwdriver
		9100-0403	S 240mm	5*130mm Pin 6 PCS 5mm Hex Wrench 5mm Screwdriver

CHANGZHOU MEDITECH TECHNOLOGY CO., LTD

ADD:9th floor, INNO LASER Building, Science & Education Town, Wujin District, Changzhou, Jiangsu, China

WhatsApp: +86-18112515727

FAX: 86-519-85855955

Email: song@orthopedic-china.com
